



# Thermoelectric Controller Peltier System



Based on Peltier effect, the thermoelectric controller is used to analyze and measure samples under the constant temperature. It is applicable for Spectro UV-VIS-2700 and UVS-2800, Spectro UVD-2950, UVD-3000, UVD-3200, UVD-3400 and UVD-3500.

Basic configuration: a temperature controller and a constant-temperature cell holder. Types of accessories:

- (1) single-stage constant-temperature device.
- (2) multi-stage constant-temperature device.

Labomed, Inc. is certified by ISO-9001-2000, has CE Conformity and is FDA Licensed.

## Technical Specifications

Temperature range:	15°C-55°C	Adjustment mode:	Continuous adjustment
Temperature accuracy:	0.1°C	Power supply:	85VAC-265VAC, 50-60Hz
Temperature resolution:	0.0625°C	Ambient temperature:	0-40°C
Display mode:	LCD digital display	Relative humidity:	<80%